

Exp. No. 5	Experiment/Subject	Date 10/10/2023
Name Emmean Kai Iash Ramush	Lab Partner	Locker/Desk No. 24
		Course & Section No. 213 203

Empirical formula of zinc iodide.

Reference: "An Experiment in Thinking of Scientifically", D.J. Sandella Chem. Ed. 69:933 (1992).

Purpose: To determine the Empirical formula of Zinc iodide, Binary Compound from product that been produced by react between Zinc and Iodide.

Materials - Electronic balance, ~~stand~~ volumic flask, Breaker, hood, zinc, Iodine, Electro heater.

Procedure.

1. Record the mass of ^{dry} (PED) directly on notebook.
2. Weight the 2.00 grams of zinc and add PED. Weight PED with 2.00 grams.
3. Calculate the different mass of PED before and after.
4. Weight Vial containing the iodine crystal (Put the cap on)
5. Transfer the iodine crystal into (PED) that have zinc
6. Weight the (PED) with zinc and iodine.
7. Determined the zinc crystal.
8. take (PED) with iodide to fume hood and observe the physical changes.
9. Add 5ml distilled water to (PED) that contained iodide and swirl

Signature	Date	Witness/TA	Date
-----------	------	------------	------

Exp. No.	Experiment/Subject	Date
Name	Lab Partner	Locker/ Desk No.
		Course & Section No.

observe it.

10. In the fume hood Add 6-M acetic solution to iodine and water
11. Observe an chemical change (it take 10 minutes)
12. ~~also~~ clean (PED) with zinc iodide three time with 1mL portion distilled water with two or three drop of 6-M acetic solution
13. dry the PED by and place it on hot plate (if dry oven is full)
14. dry zinc and rest until it reach room temperature
15. Weight (PED) with dry and unreacted zinc
16. determine the mass of unreacted zinc
17. determine the mass of zinc that react to iodide.

Signature	Date	Witness/TA	Date
-----------	------	------------	------

Exp. No.	Experiment/Subject	Date
Name	Lab Partner	Locker/ Desk No.
		Course & Section No.

Mass of Evaporating Dish (g) 41.786g
 Mass of Evaporating Dish + Zinc (g) 42.787g
 Mass of Zinc (g) 1.001g
 Mass of Vial filled with iodine (g) 6.765
 Mass of Vial Empty (g) 5.760
 Mass of iodine (g) 1.005
 Mass of Evaporating Dish + Unreacted zinc 42.522g
 Mass of unreacted zinc (g) 0.736g

Calculation

Mass of Zinc

$$41.786 - 42.787 = -1.001$$

Mass of iodine

$$6.765 - 5.760 = 1.005$$

Signature	Date	Witness/TA	Date
		ND	